



Phone: 1 800 222-4999
 Fax: 1 800 645-2519

Account Information

Account Number	Date
Account Name	Contact
Phone	Fax
Patient ID	
Prescribing Physician	

Circular Knit Compression Garment Information (without seam)

Compression: 20-30 mmHg Model: Dynamic (Varin Soft-in)
 30-40 mmHg Dynamic Silver (Varin Soft-in)


Style: C-F below knee F-G above knee H-D Pirogoff
 Silicone Border Hip Attachment

K-T Body Part (worn with F-G) Slip-on Hook and Loop Closure
 20-30 mmHg 30-40 mmHg

Extremity: Right Left Quantity: _____ Pieces

Comments:

Custom Measurement Form for Stump Shrinkers



Measure circumferences beginning 5 cm from the distal end of the stump.

S = Distal end of the stump. Measure all lengths from 'S'.